



OPHTHALMOLOGY AND OTOLARYNGOLOGY

Ophthalmology and Otolaryngology

Ultracell® Eye Products

Aspen Surgical manufactures a complete line of sponge products used for management of fluids during ophthalmic procedures. Our PVA and cellulose sponge products are used during cataract and other refractive surgeries. For LASIK procedures, our UltraSpear™, ultra clean PVA surgical eye spears, are ideal for smoothing the flap and manipulating tissue. We also offer corneal light shields, drains, and wicks in various sizes and capacities.

Aspen Surgical's Specially Formulated PVA Material:

- Ultracell® PVA is a synthetic sponge composed of polyvinyl alcohol
- Specially formulated PVA material is virtually lint-free
- Made of ultra fine pores
- Continuous open pore construction creates strength as well as softness
- Highly absorbent
- Instantaneous expansion
- Biocompatible

LASIK Products

Aspen Surgical offers products specially formulated for use during LASIK surgery. We recommend the following line of products for use during LASIK surgery.

LASIK PVA Spears and Disk

Aspen Surgical offers UltraSpear™ LASIK eye spears. The UltraSpear™ material features our exclusive patent-pending particulate extraction process to make it the cleanest eye spear we offer. This means they're virtually particulate free. UltraSpear™ Classic is a traditional compressed shape and is the fastest expanding spear available. Our pre-expanded UltraSpear™ features instantaneous wicking. The pre-expanded tip allows surgeons to adjust the firmness for better control during tissue manipulation. Our LASIK Disk is used for blotting and drying the corneal flap and bed. Latex-free and intended for single-use.



Part No.	Description	Quantity
40800	UltraSpear™ Expanded	2/Pack, 20/Box
40805	UltraSpear™ Expanded	10/Pack, 10/Box
40813	UltraSpear™ Classic	3/Pack, 20/Box
40815	UltraSpear™ Classic	10/Pack, 10/Box
40804	LASIK Disk	1/Pack, 20/Box

LASIK Drains & Shields

Aspen Surgical also offers a family of other Ultracell® brand PVA products specially designed for the fluid management needs unique to LASIK surgery. Products include drains and shields for both nasal and superior hinges. Latex-free and intended for single-use.



Part No.	Description	Quantity
40820	LASIK Shield (7mm flat spot)	1/Pack, 20/Box
40821	LASIK Shield (6mm flat spot)	1/Pack, 20/Box
40822	LASIK Shield (4mm & 8mm flat spots)	1/Pack, 20/Box
40830	LASIK Drain (11.5mm opening)	1/Pack, 20/Box
40831	LASIK Drain (10.5mm opening)	1/Pack, 20/Box

Corneal Marking Pad

Aspen's Corneal Marking Pad can be used with many refractive/corneal marking instruments. The pre-inked, ready-to-use gentian violet marking pad saves prep time and allows you to focus on your patient. Sterile, single-use, and latex-free.

Part No.	Description	Quantity
40490	Aspen Corneal Marking Pad	1/Pack, 10/Box



Cataract and Refractive Surgery Products

Aspen Surgical provides a complete line of Ultracell® brand surgical eye spears including both cellulose and PVA spears used in the rapid recovery of fluids during cataract and other refractive surgeries.

PVA & Cellulose Eye Spears

Aspen Surgical manufactures both PVA and cellulose eye spears for use in ophthalmic surgery. The cellulose spear is the standard spear used in cataract surgery. The PVA eye spear is used in both cataract and LASIK surgeries and has less particulate than the cellulose spear. Latex-free and intended for single-use.




Part No.	Description	Quantity
40400-8	PVA Spears	5/Pack, 20/Box
40405-8	PVA Spears	10/Pack, 20/Box
40410	Cellulose Spears	5/Pack, 20/Box
40413	Cellulose Spears	6/Pack, 24/Box
40415	Cellulose Spears	10/Pack, 20/Box

Part No.	Description	Quantity
40401	PVA Points	5/Pack, 20/Box
40406	PVA Points	10/Pack, 20/Box
40411	Cellulose Points	5/Pack, 20/Box
40416	Cellulose Points	10/Pack, 20/Box
40412	Cellulose Points and Strips	2/Pack, 6/Pack, 20/Box

Other Cataract and Refractive Surgery Products

Aspen Surgical provides a variety of Ultracell® PVA products including eye drains, wicks, shields, ring, and pledgets to aid in the continuous collection and control of excessive fluids during surgery. Drains include a fluid collection bag, while wicks are the PVA strip only. Latex-free and intended for single-use.

Part No.	Description	Quantity
40420	Corneal Light Shield (7mm)	1/Pack, 20/Box
40421	Corneal Light Shield (8mm)	1/Pack, 20/Box
40422	Corneal Light Shield (9mm)	1/Pack, 20/Box
40430	Eye Drain (80cc capacity)	1/Pack, 20/Box
40431	Eye Wick (4mm x 170mm)	2/Pack, 10/Box
40435	Eye Drain (400cc capacity)	1/Pack, 10/Box
40440	Bloomberg Ophthalmic Ring	1/Pack, 20/Box
40450	Eye Kit (1 Wipe/Wick, 1 Light Shield)	1/Pack, 20/Box
40451	Cataract Kit (5C. Spear, Lgt Shld, Eye Drn, Inrt Wp)	1/Pack, 20/Box
40825	Eye Pledget	10/Pack, 20/Box




Ultracell® ENT Products

Aspen Surgical manufactures a complete line of medical sponge products for otolaryngology including the nasal, sinus, epistaxis, and otology sub-specialties. Along with our products for instrument cleaning and care, all Ultracell® ENT sponge products are produced using specially formulated PVA material making them both highly absorbent and soft. For the surgeon, this means high performance. For the patient, this means increased comfort.

Nasal Packs

Aspen Surgical provides a complete line of Ultracell® brand PVA sponges for postoperative packing following nasal surgery.

- Used following septoplasty, rhinoplasty, and turbinectomy procedures
- May aid in holding septum firmly in place
- Available in a range of sizes: 3.0cm pediatric to 10.0cm
- Select sizes available in an anatomical shape
- Popular sizes also available in our patented layered design
- Excellent alternative to cumbersome gauze packing
- Strung for safety
- Latex-free
- Single-use



Regular Nasal

Part No.	Description	Quantity
10096-C	3cm Pediatric Nasal Slim Pack	1/Pack, 20/Box
10100-C	4.5cm Nasal Pack	1/Pack, 20/Box
10101-C	4.5cm Nasal Slim Pack	1/Pack, 20/Box
10105-C	6.0cm Nasal Pack	1/Pack, 20/Box
10106-C	6.0cm Nasal Slim Pack	1/Pack, 20/Box
10110-C	8.0cm Nasal Pack	1/Pack, 20/Box
10111-C	8.0cm Nasal Slim Pack	1/Pack, 20/Box
10115-C	10.0cm Nasal Pack	1/Pack, 20/Box
10120-C	8.0cm Anatomical Pack	1/Pack, 20/Box
10125-C	8.0cm Anatomical Slim Pack	1/Pack, 20/Box
10126-C	5.5cm Pediatric Anatomical Pack	1/Pack, 20/Box
10127-C	5.5cm Pediatric Anatomical Slim Pack	1/Pack, 20/Box



Layered Nasal

Part No.	Description	Quantity
10100-L	4.5cm Nasal Pack, Layered	1/Pack, 20/Box
10105-L	6.0cm Nasal Pack, Layered	1/Pack, 20/Box
10110-L	8.0cm Nasal Pack, Layered	1/Pack, 20/Box
10115-L	10.0cm Nasal Pack, Layered	1/Pack, 20/Box
10120-L	8.0cm Anatomical Pack Layered	1/Pack, 20/Box

UltraPlus 5000™ Nasal Pack

Aspen's UltraPlus 5000™ Nasal Pack features an 8.0cm nasal pack inside custom wrapping to help ease the pain of removing the nasal pack. Latex-free, intended for single-use, and strung for safety.




Part No.	Description	Quantity
50110	UltraPlus 5000 8.0cm Nasal Pack	1/Pack, 10/Box
50111	UltraPlus 5000 8.0cm Nasal Slim Pack	1/Pack, 10/Box

Sinus Packs

Aspen Surgical offers a complete line of Ultracell® brand sinus packing specifically designed for use following sinus surgeries.

- Ideal packing following FESS (Functional Endoscopic Sinus Surgery)
- Used for post-operative packing or dressing of sinus cavity
- Available in a range of lengths and thicknesses
- Patented layered design allows for easy manipulation
- Strung for safety
- Latex-free
- Single-use



Part No.	Description	Quantity
20200-C	2.5cm Sinus Pack	2/Pack, 10/Box
20200-L	2.5cm Sinus Pack, Layered	2/Pack, 10/Box
20201-C	2.5cm Sinus Pack Slim	2/Pack, 10/Box
20202-C	2.5cm Sinus Pack, Ultra Slim	2/Pack, 10/Box
20210-C	3.5cm Sinus Pack	2/Pack, 10/Box
20210-L	3.5cm Sinus Pack, Layered	2/Pack, 10/Box
20211-C	3.5cm Sinus Pack, Slim	2/Pack, 10/Box
20212-C	3.5cm Sinus Pack, Ultra Slim	2/Pack, 10/Box

Epistaxis Products

The Ultracell® brand of epistaxis sponges provides highly absorbent treatment of nosebleeds.

- Fluid absorption causes the packing to expand producing a tamponading effect
- Available in a variety of sizes: 3.0cm anterior to 10.0cm posterior
- All are strung for safety
- Ideal for use in emergency rooms, clinics, and doctor's offices
- Latex-free
- Single-use



Part No.	Description	Quantity
10140-C	10.0cm Epistaxis Pack	1/Pack, 10/Box
10140-L	10.0cm Epistaxis Pack, Layered	1/Pack, 10/Box
10141-C	10.0cm Epistaxis Slim Pack	1/Pack, 10/Box
10145-C	8.0cm Epistaxis Pack	1/Pack, 10/Box
10146-C	8.0cm Epistaxis Slim Pack	1/Pack, 10/Box
10150-C	5.5cm Epistaxis Pack	1/Pack, 10/Box
10150-L	5.5cm Epistaxis Pack, Layered	1/Pack, 10/Box
10151-C	5.5cm Epistaxis Slim Pack	1/Pack, 10/Box
10160-C	3.0cm Anterior Epistaxis Pack	1/Pack, 10/Box

Ophthalmology and Otolaryngology

Otology Products

Aspen Surgical provides both ear wicks and ear packs.

- Much smaller, compressed form allows easy insertion into swollen ear canal with minimal patient discomfort
- Pediatric sizes available



Ear Wicks

- Perfect for use with treatments for otitis externa
- Help maintain antibiotics in the ear canal
- Also available fenestrated
- Latex-free
- Single-use


Part No.	Description	Quantity
30300-C	Ear Wick (9mm x 15mm)	1/Pack, 50/Box
30301-C	Pediatric Ear Wick (7mm x 12mm)	1/Pack, 50/Box
30305-C	Ear Wick, Fenestrated (9mm x 15mm)	1/Pack, 50/Box
30306-C	Pediatric Ear Wick, Fenestrated (7mm x 12mm)	1/Pack, 50/Box
30310-C	Ear Wick, Long (9mm x 24mm)	1/Pack, 10/Box
30311-C	Pediatric Ear Wick, Long (7mm x 20mm)	1/Pack, 10/Box
30315-C	Ear Wick, Fenestrated, Long (9mm x 24mm)	1/Pack, 10/Box
30316-C	Pediatric Ear Wick, Fenestrated, Long (7mm x 20mm)	1/Pack, 10/Box



Ear Packs

- Used following otic surgery
- Provides post-surgical packing and dressing of ear canal
- Latex-free
- Single-use

Part No.	Description	Quantity
30320-C	Ear Pack (9mm x 15mm)	1/Pack, 10/Box
30330-C	Ear Pack Large (12mm x 15mm)	1/Pack, 10/Box
30340-C	Ear Pack Large, Long (12mm x 24mm) No string	1/Pack, 10/Box



Micro Spheres

For tissue manipulation during ear surgery. Latex-free and intended for single-use.

Part No.	Description	Quantity
30380-C	3mm Micro Sphere	10/Pack, 20/Box
30381-C	4mm Micro Sphere	10/Pack, 20/Box
30382-C	5mm Micro Sphere	10/Pack, 20/Box

Specialty Products

Aspen Surgical also provides specialty PVA products. Included are PVA wipes used in the care and cleaning of surgical instruments. We also manufacture a PVA surgical headband. Wipe and cleaning block intended for single-use. Headband can be washed and re-used. All are latex-free.




Part No.	Description	Quantity
40462	Diamond Knife Cleaning Block	1/Pack, 20/Box
90900	Instrument Wipe (2mm)(3"x3"x2mm)	1/Pack, 20/Box
90910	Instrument Wipe (1mm)	1/Pack, 20/Box
90950	Headband	1/Pack, 20/Box

Ophthalmic Needles

Aspen Surgical offers three types of needles for Ophthalmic procedures. The reverse cutting needles are ideal for ocuplastic surgery as the needle allows easy passage through the epidermis. The taper point produces very small holes and is useful in iris repair. Richard-Allan surgical needles have unsurpassed strength, sharpness, and smoothness for effortless penetration and reduced tissue trauma. Each needle is offered in several sizes as listed below, and each package contains 2 sterile needles, 72 packages per box. Latex-free and intended for single-use.

2123

Eye Needle ½ Circle Reverse Cutting



212301 0.024 in
1.024 in



212303 0.021 in
0.748 in



212302 0.024 in
0.866 in




212304 0.021 in
0.630 in



2124

Eye Needle ¾ Circle Reverse Cutting



212401 0.024 in
1.260 in



212404 0.021 in
0.866 in



212402 0.024 in
1.063 in



212405 0.021 in
0.748 in



212403 0.024 in
1.063 in




212406 0.021 in
0.630 in



2125

Eye Needle ½ Circle Taper Point



212501 0.024 in
1.063 in



212503 0.021 in
0.748 in



212502 0.024 in
0.866 in



212504 0.021 in
0.630 in



In addition to Ultracell®, Aspen Surgical carries a full line of O.R. Accessory, Woundcare, Instrument Care, and Patient and Staff Safety products. Please call us to request one of the catalogs above, or to share your new product ideas!



6945 Southbelt Dr SE
Caledonia, MI 49316 USA
Phone: 616.698.7100
888.364.7004
Fax: 616.698.9281
Web: aspensurgical.com

International Sales:
NETWork Medical Products Ltd.
Phoenix Business Centre, Ripon, UK
Tel: +44 (0)1765 609555
Fax: +44 (0)1765 608476
www.networkmedical.co.uk