



Duckworth & Kent Ltd.
titanium surgical instrument manufacturer

at the leading edge

NEWSLETTER

MAY 2018



CORE VALUES
Collection

E-RANGE
Instruments

Essential
CATARACT SET

Find us on



E-RANGE Instruments




Innovation is at the heart of our philosophy and is driven by our enduring commitment to provide surgeons with only the highest quality reusable titanium ophthalmic instrumentation.

We are proud to show our latest line of instrumentation, the "E-Range".

This range represents the application of cutting edge technology fused with over 50 years of engineering knowledge. These instruments offer surgeons an elegant, well balanced and ergonomically designed handle whilst still providing our high quality, consistent and durable functionality.

2-716GE-1 & 2-716GE-1S Inamura Flat Handle Capsulorhexis Forceps 1.8mm Incision




2-716GE-1



2-716GE-1S


- Sharp pointed serrated tips
- Cross action tips, 1.5mm width at pivot box
- Curved shaft, tip to pivot point **10.5mm**
- Marks on shaft 2.5mm and 5mm denote desired size of capsulorhexis
- Tips angled 45° from shaft
- Flat handle, length 111mm = 2-716GE-1
- Short flat handle, length 86mm = 2-716GE-1S

2-716GER8 Inamura Round Handle Capsulorhexis Forceps 1.8mm Incision



- Sharp pointed serrated tips
- Cross action tips, 1.5mm width at pivot box
- Curved shaft, tip to pivot point **9.5mm**
- Marks on shaft 2.5mm and 5mm denote desired size of capsulorhexis
- Tips angled 45° from shaft
- Round handle, length 110mm

2-716GER8-2 Inamura Round Handle Capsulorhexis Forceps 1.8mm Incision



- Sharp pointed serrated tips
- Cross action tips, 1.5mm width at pivot box
- Curved shaft, tip to pivot point **8.5mm**
- Marks on shaft 2.5mm and 5mm denote desired size of capsulorhexis
- Tips angled 45° from shaft
- Round handle, length 110mm

Visit www.duckworth-and-kent.com to view our full range of E-Range instruments.

Thank you

Duckworth & Kent would like to thank all those who visited our distributors booth at the ASCRS in Washington DC. We hope that you enjoyed your visit.


Please visit our events page at www.duckworth-and-kent.com/events to keep up to date with all future meetings.



ICL Instruments Overview

2-770N

Deitz ICL Loading Forceps




- Inside jaws textured, 4mm tip angle
- 20° angled shafts, tip to angle length 14mm
- Highly polished outer jaws
- Flat handle, length 83mm

Extra delicate, long jaws for inserting ICL into barrel of injector cartridge.



2-896-3

Small Incision ICL Manipulating Forceps, Angled Flat Disc Tips



- Flat tip, 1.9mm wide disc with platform
- Upper disc coloured gold for orientation reference
- Tips angled at 15°
- Corresponding gold star on handle, relates to gold tip



- Squeeze action handle activates both jaws
- 20 gauge, 32.5mm long straight shaft
- Round squeeze handle, length 139mm

Visit www.duckworth-and-kent.com to view our full range of ICL instruments.


Premium Quality at Competitive Prices

As a continuation of the process to utilise innovative engineering, Duckworth & Kent is very pleased to present to you its "Core Values Collection"

This diverse range of reusable instruments combines our renowned precision and quality at an exceptional price.

Scan the QR code to view our dedicated "Core Values Collection" leaflet.

CORE VALUES
Collection




Visit www.duckworth-and-kent.com to view our full range of Core Value instruments.

Essential Cataract Set Selection

8-652S

Short Handle Irrigation Handpiece




- 23 gauge, 0.65mm tube diameter
- Two 0.4mm side irrigation ports

- Curved shaft, tube length 15mm
- Round handle, length 80mm



8-657S

Short Handle Aspiration Handpiece



- 23 gauge, 0.65mm tube diameter, textured tip
- 0.3mm aspiration port

- Curved shaft, tube length 11mm
- Round handle, length 77mm



View our complete set of Essential Cataract instruments →

Essential CATARACT SET

Visit www.duckworth-and-kent.com to view our full range including our cost calculator.

The decision to choose our reusable instruments over single-use alternatives is a simple one. Not only are you investing in world-renowned high quality titanium alloy products, but also in instruments that have less impact on the environment and are more cost-beneficial in the long term.

When you invest in a Duckworth & Kent instrument you are purchasing a product that will serve you impeccably for innumerable years and countless procedures.

When compared with single-use alternatives, Duckworth & Kent instruments are not only more cost-beneficial long term, but have a lower impact on the environment. See how much you can save by investing in our innovative precision engineered reusable titanium product range.



Cost Calculator for Essential Cataract Set



If you are interested in any of our product range or have a general enquiry please contact anyone of our dedicated customer service team members who are ready to help.



Maria
Customer Services
Supervisor



Kelly
Customer Services



mariagente@
kellyobrien@
nickyguthrie@ } duckworth-and-kent.co.uk

Find us on



© May 2018 Duckworth & Kent

Email: info@duckworth-and-kent.co.uk

Tel: +44 (0)1462 893254

Fax: +44 (0)1462 896288

www.duckworth-and-kent.com

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.