



Advanced. Versatile. Different.

Stellaris PC
Product Catalogue


BAUSCH + LOMB
Stellaris[®] PC
Vision Enhancement System

Stellaris PC Hardware + Accessories

BL1433
Stellaris PC Combined System
with 1 Xenon + 1 Xenon-Mercury
Lamp Assembly  = 1



BL2393
Stellaris PC Foot Pedal  = 1



- BL4390 Foot Control Battery
- BL4391 Foot Control Wall Charger
- BL4393 Foot Control Charging Cradle
- BL4394 Foot Control Power Cable

Stellaris PC Hardware + Accessories

BL3270
Fragmentation Handpiece  = 1



BL3170
Stellaris Phaco Handpiece  = 1



BL3320
Xenon lamp Replacement  = 1



BL3321
Mercury Lamp Replacement  = 1



BL5275
20GA Fragmentation Needle
15 deg single use  = 6



BL4363
Stellaris PC Zero Level Bottle
Hanger - Dual Bottle  = 1



BL2300
Stellaris PC/Stellaris Multi-
Media Center (MMC)  = 1



BL3253EN
Stellaris PC User Manual, English  = 1

- BL3253DE German
- BL3253FR French
- BL3253ES Spanish
- BL3253IT Italian
- BL3253NL Dutch
- BL3253PT Portuguese
- BL3253SV Swedish
- BL3253NO Norwegian



BL4360
Stellaris PC/Stellaris Mayo Tray Adapter

Please contact your local representative to
order a Mayo Tray.

Contact details can be found on the
back page.




Stellaris PC Procedural Packs Combined with AFI (Air Forced Infusion)

BL5520
20G Combined Pack




- Pack Includes:**
- A 20G priming cup subassembly
 - B Scleral plug package - 19 & 20
 - C 20G light pipe
 - D 20G MVR knife



BL5523
23G Combined Pack




- Pack Includes:**
- A 23G priming cup subassembly
 - B 23G light pipe
 - C 23G ESA System
 - D 23G Scleral plugs



BL5525
25G Combined Pack




- Pack Includes:**
- A 25G priming cup subassembly
 - B 25G light pipe
 - C 25G ESA System
 - D 25G Scleral plugs



All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 10 cc syringe
- 6 Cassette stopper
- 7 Silicone test chamber
- 8 Vitrectomy Cutter Handle Extension
- 9 Male luer coupler
- 10 Irrigation sleeves x 2
- 11 Needle wrench
- 12 Pre-piercing spike assembly




Stellaris PC Procedural Packs Combined

BL5320
20G Combined Pack




- Pack Includes:**
- A 20G priming cup subassembly
 - B Scleral plug package - 19 & 20
 - C 20G light pipe
 - D 20G MVR knife



BL5323
23G Combined Pack




- Pack Includes:**
- A 23G priming cup subassembly
 - B 23G light pipe
 - C 23G ESA System
 - D 23G Scleral plugs



BL5325
25G Combined Pack




- Pack Includes:**
- A 25G priming cup subassembly
 - B 25G light pipe
 - C 25G ESA System
 - D 25G Scleral plugs



All the above packs come with the items below.

- 1 Cassette with combined gravity tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 10 cc syringe
- 6 Cassette stopper
- 7 Silicone test chamber
- 8 Vitrectomy Cutter Handle Extension
- 9 Male luer coupler
- 10 Irrigation sleeves x 2
- 11 Needle wrench




Stellaris PC Procedural Packs

Posterior with AFI (Air Forced Infusion)

BL5420
20G Posterior Pack




- Pack Includes:
- A 20G priming cup subassembly
 - B Scleral plug package - 19 & 20
 - C 20G light pipe
 - D 20G MVR knife



BL5423
23G Posterior Pack




- Pack Includes:
- A 23G priming cup subassembly
 - B 23G light pipe
 - C 23G ESA System
 - D 23G Scleral plugs



BL5425
25G Posterior Pack




- Pack Includes:
- A 25G priming cup subassembly
 - B 25G light pipe
 - C 25G ESA System
 - D 25G Scleral plugs



All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 10 cc syringe
- 6 Cassette stopper
- 7 Silicone test chamber
- 8 Vitrectomy Cutter Handle Extension
- 9 Male luer coupler
- 10 Pre-piercing spike assembly




Stellaris PC Procedural Packs

Posterior

BL5220
20G Posterior Pack




- Pack Includes:
- A 20G priming cup subassembly
 - B Scleral plug package - 19 & 20
 - C 20G light pipe
 - D 20G MVR knife



BL5223
23G Posterior Pack




- Pack Includes:
- A 23G priming cup subassembly
 - B 23G light pipe
 - C 23G ESA System
 - D 23G Scleral plugs



BL5225
25G Posterior Pack




- Pack Includes:
- A 25G priming cup subassembly
 - B 25G light pipe
 - C 25G ESA System
 - D 25G Scleral plugs



All the above packs come with the items below.

- 1 Cassette with gravity tube assembly
- 2 Second Aspiration Line
- 3 Irrigation Line
- 4 Drapes for Mayo tray, screen and remote control
- 5 10 cc syringe
- 6 Cassette stopper
- 7 Silicone test chamber
- 8 Vitrectomy Cutter Handle Extension
- 9 Male luer coupler



VRL100
DK-Line® 5ml



VRL110
DK-Line® 7ml



VRL200
Okta-Line® 5ml



VRL500
Oxane® 1300, 10ml
Syringe



VRL600
Oxane® 5700, 10ml
Syringe



VRL700
Oxane® HD, 10ml Syringe



CX5710
Viscous Fluid Extraction
Kits



CX5720
Viscous Fluid Injection Kits



Forceps


65023S
23G Single Use
Endgripping VR Forceps




65024S
23G Single Use
Endgripping VR Forceps



65025S
25G Single Use
Endgripping Forceps



65026S
25G Single Use
Endgripping Forceps



SUI03
23G Asymmetric Peeling,
ILM, Forceps




SUI04
23G Asymmetric Peeling,
ILM, Forceps



SUI05
25G Asymmetric Peeling,
ILM, Forceps



SUI06
25G Asymmetric Peeling,
ILM, Forceps



Membrane Picks

SUK01
23G Membrane Pick



SUK02
25G Membrane Pick



Bipolar Eraser

E7925 25G Bipolar Eraser, Blunt, 2-pin connector	= 12 STERILE		
E7924 25G Bipolar Eraser, Tapered, 2-pin connector	= 12 STERILE		
E7931 23G Bipolar Eraser, Tapered, 2-pin connector	= 12 STERILE		
E7928 Scleral Bipolar Eraser, Lemosa type connector	= 12 STERILE		
E7929 20G Bipolar Eraser, Tapered, Lemosa type connector	= 12 STERILE		
E7930 23G Bipolar Eraser, Tapered, Lemosa type connector	= 12 STERILE		
E7926 25G Bipolar Eraser, Tapered, Lemosa type connector	= 12 STERILE		
E7927 25G Bipolar Eraser, Blunt, Lemosa type connector	= 2 STERILE		

Aspiration Handpieces

SUJ01 Passive Aspiration Handpiece with 23G Soft Tip cannula.	= 12 STERILE		
SUJ04 Passive Aspiration Handpiece with 23G Soft Tip cannula.	= 2 STERILE		
SUJ02 Passive Aspiration Handpiece with 25G Softtip cannula.	= 12 STERILE		
SUJ05 Passive Aspiration Handpiece with 25G Softtip cannula.	= 2 STERILE		
BL7823 23G Soft Tip Cannula	= 10 STERILE		
BL7725 25G Soft Tip Cannula	= 10 STERILE		

Knives


E7519 19G MVR Blade	= 6 STERILE		A: 31.8mm B: 1.5mm C: 3.1mm
E7520 20G MVR Blade	= 6 STERILE		A: 31.8mm B: 1.2mm C: 2.4mm
E7525 25G MVR Blade	= 6 STERILE		A: 31.8mm B: 0.5mm C: 1.0mm

Stellaris Procedure Packs (Anterior Phaco)

BL5110
Stellaris Premium Vacuum
Phacoemulsification Pack




- Pack Includes:**
- A VFM cassette assembly and I/A set
 - B Auxillary drape
 - C Screen drape
 - D Remote control drape
 - E Irrigation sleeves x 2
 - F Plug
 - G Test chamber
 - H Needle wrench



BL5112
Stellaris Optimised Stability
Vacuum Pack




- Pack Includes:**
- A VFM cassette assembly and I/A set
 - B Auxillary drape
 - C Screen drape
 - D Remote control drape
 - E Irrigation sleeves x 2
 - F Plug
 - G Test chamber
 - H Needle wrench



BL5113
Stellaris StableChamber™ Phaco
Pack with MICS needle




- Pack Includes:**
- A VFM cassette assembly and I/A set
 - B Auxillary drape
 - C Screen drape
 - D Remote control drape
 - E Irrigation sleeves x 2
 - F Plug
 - G Test chamber
 - H Needle wrench
 - I StableChamber
 - J MICS straight needle



BL5113A
Stellaris StableChamber™ Phaco
Pack with MICS angled needle



- Pack Includes:**
- A VFM cassette assembly and I/A set
 - B Auxillary drape
 - C Screen drape
 - D Remote control drape
 - E Irrigation sleeves x 2
 - F Plug
 - G Test chamber
 - H Needle wrench
 - I StableChamber
 - J MICS angled needle




Stellaris Procedure Packs (Anterior Phaco)

BL5114
Stellaris Premium Vacuum Phaco
Pack with MICS needle



- Pack Includes:**
- A VFM cassette assembly and I/A set
 - B Auxillary drape
 - C Screen drape
 - D Remote control drape
 - E Irrigation sleeves x 2
 - F Plug
 - G Test chamber
 - H Needle wrench
 - I MICS straight needle



BL5114A
Stellaris Premium Vacuum Phaco
Pack with MICS angled needle



- Pack Includes:**
- A VFM cassette assembly and I/A set
 - B Auxillary drape
 - C Screen drape
 - D Remote control drape
 - E Irrigation sleeves x 2
 - F Plug
 - G Test chamber
 - H Needle wrench
 - I MICS angled needle



BL5612
20GA High Speed
(30-5000cpm) Anterior
Vitrectomy Cutter Pack



Single Use Cataract Instruments

Manipulators

SUF01 Koch Stop & Chop Manipulator	= 12 STERILE		
SUF02 Sinsky Hook	= 12 STERILE		
SUF03 Phaco Spatula	= 12 STERILE		

Cannulas

82020S Corneal Irrigator	= 60 STERILE		
-----------------------------	-----------------	--	--

Speculum

H01 Barraquer Eye Speculum available as Per Procedure Tray component part	= 1 STERILE	
52001S Eye Speculum, adjustable	= 3 STERILE	
52000S Eye Speculum with aspiration, adjustable	= 3 STERILE	

Single Use Cataract Instruments

Forceps

SUE01 Capsulorhexis Forceps	= 8 STERILE		
SUE02 Bonn Forceps	= 8 STERILE		
SUE03 Tying Forceps	= 8 STERILE		
SUE04 McPherson Tying Forceps, angled	= 8 STERILE		
SUE05 Pierse Tying Forceps	= 8 STERILE		
SUE06 MICS™ Capsulorhexis Forceps, fine angled tips, angled shafts for 1.8mm incision.	= 8 STERILE		

Knives

E7595 14-16mm clear corneal	= 6 STERILE		A: 12.6mm B: 1.6mm C: 1.7mm
E7599 15-17mm trapezoidal	= 6 STERILE		A: 14.4mm B: N/A C: 2.0mm D: 2.0mm E: N/A
E7596 15-2.0mm trapezoidal	= 6 STERILE		A: 14.4mm B: 2.0mm C: 1.8mm D: 1.5mm E: N/A
E7600 16-18mm trapezoidal	= 6 STERILE		A: 14.4mm B: 1.8mm C: 1.8mm D: 1.6mm E: 2.95mm

Single Use Cataract Instruments

E7551A 2.65mm, angled up	= 6 STERILE		A: 14.6mm B: 2.6mm C: 3.9mm
E7559A 2.75mm, angled up	= 6 STERILE		A: 14.4mm B: 2.75mm C: 3.2mm
E7553A 2.85mm, angled up	= 6 STERILE		A: 14.1mm B: 2.8mm C: 3.3mm
E7515 15° sideport	= 6 STERILE		

Stellaris® System Coaxial I/A Handpieces

85794ST I/A handpiece curved shaft tip, with sleeve for 1.8mm incision	= 12 STERILE		
85795ST I/A handpiece 45° angled tip, with sleeve for 1.8mm incision	= 12 STERILE		
85910ST Silicone I/A handpiece 45° incisions: 2.2mm-2.8mm	= 12 STERILE		
85783ST 17G (Ø 1.5mm) I/A Handpiece, curved shaft	= 12 STERILE		
85784ST 17G (Ø 1.5mm) I/A Handpiece straight shaft 45° angled tip 0.3mm aspiration port	= 12 STERILE		
85786ST 16G (Ø 1.65mm) I/A Handpiece straight shaft, 90° angled tip, 0.3mm aspiration port	= 12 STERILE		

Single Use Cataract Instruments

Stellaris® System Biaxial I/A Handpieces

85902ST19 19G (Ø 1.07mm) Irrigation handpiece, curved shaft, smooth tip, 2 lateral irrigation port of 0.6mm	= 12 STERILE		
85901S19 19G (Ø 1.07mm) Aspiration handpiece, curved shaft, roughened tip, 0.3mm aspiration port	= 12 STERILE		
85787ST 21G (Ø 0.8mm) Irrigation handpiece, curved shaft, smooth tip, 2 lateral irrigation port of 0.6mm	= 12 STERILE		
85780S 21G (Ø 0.8mm) Aspiration handpiece, curved shaft, roughened tip, 0.3mm aspiration port	= 12 STERILE		

Stellaris® System Irrigating Chopper


85903ST19 19G (Ø 1.07mm) Nagahara type tip, end out	= 12 STERILE		
--	-----------------	--	--

Coaxial I/A Handpieces

85910S Silicone I/A handpiece 45° tip, incisions: 2.2mm-2.8mm	= 12 STERILE		
85782S 17G (Ø 1.5mm) I/A handpiece, straight shaft and tip	= 12 STERILE		
85783S 17G (Ø 1.5mm) I/A handpiece, curved shaft	= 12 STERILE		
85784S 17G (Ø 1.5mm) I/A handpiece, straight shaft, tip 45°	= 12 STERILE		

Single Use Cataract Instruments

85785S
I/A handpiece straight shaft with sleeve for standard coaxial surgical technique. 20G (0.9mm) aspiration cannula



85786S
16G (Ø 1.65mm) I/A handpiece straight shaft, tip 90°



Biaxial I/A Handpieces Aspiration

85901S19
Aspiration handpiece 19G (Ø 1.07mm)



85901S20
Aspiration handpiece 20G (Ø 0.89mm)



curved shaft roughened tip

85780S
Aspiration handpiece 21G (Ø 0.8mm)



85790S
Aspiration handpiece 22G (Ø 0.7mm)



85792S
Aspiration handpiece 23G (Ø 0.6mm)



Biaxial I/A Handpieces Irrigation

85902S19
19G (Ø 1.07mm) Irrigation handpiece



85902S20
20G Irrigation handpiece (Ø 0.89mm)



curved shaft smooth tip


85787S
21G Irrigation handpiece (Ø 0.8mm)



Single Use Cataract Instruments

Biaxial I/A Handpieces Irrigation


85781S
21G Irrigation handpiece (Ø 0.8mm) roughened tip



85788S
21G Irrigation handpiece (Ø 0.8mm) central port



85789S
22G Irrigation handpiece (Ø 0.7mm) central port



85791S
23G Irrigation handpiece (Ø 0.6mm)



Irrigating Choppers for Biaxial MICS™ (B-MICS) 19&20G

85903S19
19G (Ø 1.07mm) Nagahara type tip, end out



85903S20
20G (Ø 0.89mm) Nagahara type tip, end out



TSV Platform

The TSV Platform
The exclusive minimally invasive solution for
Vitreotomy and Combined Procedures.



Equipment: Stellaris PC

The new Stellaris PC System now enables high performance 5K sutureless vitrectomy with optimum surgeon control and simple functionality for every procedure.



Accessories

Combined and Posterior packs contain all disposable components, from the vitreous cutter to the entry site alignment system (ESA), needed for 20G, 23G and 25G surgery.



Instruments

The Storz® Ophthalmics brand offers a wide range of single use instruments required to perform key vitreoretinal procedures.



Support

Peace of mind for the OR team with a 24 service support package and fully certified sales personnel



Endotamponades

Vitreoretinal tamponades, Dk-line®, Okta-line®, Oxane® and Oxane Hd® ranges help in providing superior clinical outcomes in the treatment of retinal disorders.

18 MICS™

The MICS™ 1.8 Platform
The exclusive minimally Invasive solution for Cataract
and Combined Procedures



Stellaris® MICS

The Stellaris Vision Enhancement System delivers the critical capabilities required: for 1.8 mm MICS: superior fluidics control, efficient cutting dynamics and the ergonomic insights that make the system easy to use.



Akreos® MICS

The Akreos MICS Lens combines the vision benefits of an aspheric acrylic lens with the material and design characteristics that allow it to be inserted into a true 1.8 mm incision.



MICS Instruments

The MICS Platform includes precision performance instruments from Storz® Ophthalmics for both C-MICS and B-MICS procedures.
Find the complete range at www.storzeye.eu



MICS Viscoelastic

Amvisc® PLUS is the versatile viscoelastic that is ideal for every step of your 1.8 mm MICS procedure.



MICS Support

When you move to the MICS Platform you benefit from a comprehensive support program designed to assist cataract surgeons in developing and perfecting their MICS techniques.

Bausch + Lomb
European Office
106 London Road
Kingston-Upon-Thames
Surrey KT2 6TN, UK
T +44 20 8781 5500
F +44 20 8781 2997

Belgium/Luxembourg
T +32 3 280 82 40
F +32 3 280 82 59

Emerging Markets
T +33 4 67 12 30 30
F +33 4 67 12 30 32

France
T +33 4 67 12 30 30
F +33 4 67 12 30 31

Germany/Switzerland/Austria
T +49 30 330 93 5702
F +49 30 330 93 5712

Italy
T +39 029 148 3851

Netherlands
T +31 20 655 4555
F +31 20 655 4640

Nordic Countries
T +46 8 616 95 00
F +46 8 669 86 23

Portugal
T +35 802 203 178
F +35 808 203 179

Spain
T +34 902 381 010
F +34 902 250 310

South Africa
T +27 11 259 2600
F +27 11 259 2650

United Kingdom
T +44 20 8781 0000
F +44 20 8781 0001