

# PROTECTION

## (GOWNS / APRONS / TOWELS)



### Surgical Gown

Basic Material: Spun-lace, SMS, SMMS, SMMMS  
Color: Green, blue  
Reinforcement: None or reinforced in body and sleeves  
Packing: Individually packed in pouch  
Size: S, M, L, XL, XXL



### Coverall

Basic Material: Non-woven  
Weight: 25 ~ 70 g/m<sup>2</sup>  
Color: White, blue, green  
Design: Elastic in the hood, wrists and ankles; Close in front with zipper / Size: S, M, L, XL, XXL  
Packing: Plastic bag



### Isolation Gown

Basic Material: Non-woven  
Color: Yellow, blue, green  
Wristband: With cuff or elastic tie  
Packing : Plastic bag  
Size: S, M, L, XL, XXL



### Examination Table Cover

Material: 100% wooden or paper film laminated  
Weight: 36 g/m<sup>2</sup>  
Size: 49cm x 50m, 59cm x 50m, 50cm x 46m, 50cm x 100m  
Packing: Poly-bag and transport carton



### Disposable Aprons

Material: PE film, or non-woven laminated with PE  
Color: White  
Basic size: 80 x 125 cm, 80 x 140 cm, 80 x 160 cm  
Packing: 50 pcs/pack, or 100 pcs/poly-bag



### Disposable Washing Gloves

Material: Non-woven fabric  
Color: White, blue, green  
Packing: 50 pcs/poly-bag, bulk  
Size: 15 x 23 cm



### Hand towel

Material: 100% wooden or paper film laminated  
Weight: 18 ~ 55 g/m<sup>2</sup>  
Size: 38 x 40 cm, 30 x 30 cm, 40 x 40 cm  
Packing: Individual peel-open pack or poly-bag