

WOUND CARE DRESSING



Paraffin Gauze Dressing

Material:	100% cotton cloth, paraffin oil
Paraffin load:	Light loading,normal loading,heavy loading
Medication:	With or without chlorhexidine
Sterility:	Sterilized by gamma
Packing:	Aluminum lined pouch/plastic box/metal tin
Size:	5 x 5, 7.5 x 7.5, 10 x 10,10 x 15,10 x 40, 10 x 100, 10 x 700cm



Tracheotomy Compress

Material:	Non adherent net,absorbent fleece,water repellent viscose non-woven
Approx. Weight:	180,200,220,270,400 g/m2
Presentation:	With an expandable hole to accommodate tubes of diameter 12 mm to 16 mm
Packing:	Individual peel-open pack
Size:	6 x 7 cm, 8 x 9 cm, 9 x 10 cm, 9 x 15 cm



Gauze Swabs

Material:	100% cotton
Yarn:	40's, 32's, 21's
Mesh:	15 x 11, 20 x 12, 19 x 15,24 x 20
Size:	5 x 5cm, 7.5 x 7.5cm, 10 x 10cm,10 x 20cm
Sterility:	Non-sterile, sterilized(by gamma,ETO) packed ready for sterilization
Packing:	Peel-open pouch(single or double)/poly- bag/paper bag
X-ray:	With or without X-ray detectable thread(blue or black) woven-in or heat welded



Non-woven Swabs

Material: Viscose and polyester / **Weight:** 30,40,50 g/m²
Size: 5 x 5cm, 7.5 x 7.5cm, 10 x 10cm,10 x 20 cm
Ply: 4ply, 6ply, 8ply
Sterility: Non-sterile, sterilized(by gamma,ETO) packed ready for sterilization
Packing: Peel-open pouch/paper bag



Gauze Balls

Material: 100% cotton / **Yarn:** 40's, 32's
Mesh: 24 x 20, 28 x 24, 30 x 20, 30 x 28
Diameter: 2cm,3cm,4cm,5cm,6cm
Sterility: Non-sterile, sterilized(by gamma,ETO) packed ready for sterilization
Packing: Peel-open pouch/poly- bag



Non-woven Balls

Material: Viscose and polyester
Weight: 30,40,50 g/m²
Diameter: 2cm,3cm,4cm,5cm,6cm
Sterility: Non-sterile, sterilized(by gamma,ETO) packed ready for sterilization
Packing: Peel-open pouch/poly- bag



Laparotomy Sponges

Material:	100% cotton
Yarn:	40's, 32's, 21's
Mesh:	15 x 11, 20 x 12, 19 x 15, 24 x 20, 28 x 24, 30 x 20, 30 x 28
Sterility:	Non-sterile, sterilized (by gamma, ETO) packed ready for sterilization
Packing:	Paper bag/plastic bag/single or double peel open pouch
X-ray:	With or without X-ray detectable thread (blue or black) woven-in or heat secured
Size:	20cm x 20cm, 22.5cm x 22.5cm, 30cm x 30cm, 40cm x 40cm, 45cm x 45cm
Ply:	4ply, 6ply, 8ply, 12ply



Absorbent Pad

Material:	Cellulose pulp and non-woven
Color:	White, blue or green
Sterility:	Non-sterile, or sterilized by gamma, ETO
Size:	10 x 10 cm, 10 x 15 cm, 10 x 20 cm, 10 x 30 cm 15 x 25 cm, 20 x 20 cm, 20 x 30 cm, 20 x 40 cm
Packing:	1 pc/peel-open pouch, 10 pcs/bag, bulk




Adhesive Wound Dressing

Material:	Non-woven cloth
Adhesive:	Coated polymeric adhesive
Pad:	Low adherent dressing pad
Packing:	Individually packed
Size:	5 x 7.5 cm, 8 x 10 cm, 10 x 15 cm, 10 x 20 cm, 10 x 25 cm



Low adherent dressing

Material:	100% rayon and porous PE film
Features:	Low adherence to the wound, avoiding secondary injuries upon removal
Sterility:	Gamma sterile
Size:	5 x 5 cm, 7.5 x 7.5 cm, 10 x 10 cm
Packing:	Individually packaged in pouch



Waterproof Film Dressing

Material:	Polyurethane
Adhesive:	Coated acrylic adhesive
Color:	Transparent
Packing:	Individually packed in peel-open pouch
Size:	6.0 x 7.0 cm, 10 x 12 cm, 15 x 20 cm, 12 x 25 cm



Plasters

Material:	Cotton cloth, non-woven cloth, P.E. film, PVC film
Adhesive:	Coated polymeric adhesive
Color:	White, skin color
Packing:	In sheets or individually packed
Size:	1.25cm x 12cm, 1.9cm x 5.6cm, 1.9cm x 7.2cm,
Butterfly:	4.0cm x 6.5cm



Antiseptic Wipes

Material:	Tissue impregnated with Cetrimide, Chlorhexidine
Size:	14 cm x 21 cm / 10 cm x 20 cm
Color:	White
Sterility:	Non-sterile, sterilized by gamma, ETO
Packing:	Individually packed in pouch



Alcohol Pads


Material: Non-woven impregnated with 75% Ethanol or 70% Isopropyl Alcohol

Size: 3 cm x 3 cm x 2ply

Color: White

Sterility: Non-sterile,sterilized by gamma,ETO

Packing: Individually packed in pouch



Medical Applicators

Material: Polystyrene stick and foam head

Application: Clean body area that will be subject to surgery

Properties: 155 x 3 mm plastic applicator and 40 x 40 x 25 mm foam
180 x 3 mm plastic applicator and 50 x 50 x 25 mm foam

Packing: Individual peel-open pack or polybag